

# ADDITIONAL SUPPLIES

## UNDERCUT ENHANCER PLIERS

### CRIMPERS

**\$125.00**

(each)

- Improve retention
- Adjust undercut areas & embrasures of PMT and Nylon appliances


## DIAMOND BITE STICK

### A RELIABLE AND EASY APPROACH TO FINDING AN OSA/TMJ BITE

**\$15.00**

(10 pack)

- Simple, effective and globally recognized
- Peer Reviewed



## DIAMOND BEST BITE

### HYDROPHILIC BITE REGISTRATION

**\$15.00**

(Cartridge w/Tips)

- Increased dimensional stability that holds form
- 15 sec. working time; 45sec intra-oral set time



## AQUALIZER

The Aqualizer is the most accurate and effective differential diagnostic tool available. The aqualizer allows the body to find instant functional balance by unraveling the malocclusion it is locked in. The Aqualizer is an invaluable tool to anyone who is suffering from the effects of TMD. Dental devices made from this bite are immensely effective, predictable and require remarkably few adjustments, saving valuable chair time.

**\$15.00**

(each)



## RHINOMED MUTE

### NASAL DILATOR

**\$14.99**

Trial Pack  
(1 Small, 1 Medium, 1 Large)

- Mute is a nasal dilator designed to reduce snoring by gently opening nasal passages and preventing airway collapse.
- Proven to increase airflow through the Nose by an average of 38% when compared to nasal strips
- Easy to Use, flexible, lightweight and reusable up to ten times
- Available in three sizes, each with individually adjustable nostril stents.

**\$21.99**

(3 pack - Must specify Small, Medium or Large)



## DR. B's CLEANADENT CRYSTALS

### TMJ & SLEEP ORTHOTIC & PRINTED APPLIANCES

**\$10.95**

(12 oz. Bottle)

- Cleanadent™ will destroy pathogens usually found on prosthetics
- Safe for dentures soft relines, partials, orthodontic appliances, mouth guards, retainers and all TMD & Sleep orthotics



## NOVADENT ORIGINAL & IP FORMULA

### TMJ & SLEEP ORTHOTIC & PRINTED APPLIANCES

**\$7.50**

(8 pack)

- Cleanser for all dentures, splints, retainers, and all TMJ & Sleep orthotics
- IP formula recommended for printed orthotics



## APPLIANCE SAMPLES & SUPPLIES



DIAMOND ORTHOTIC  
LABORATORY  
INTERNATIONAL

# OSA DEVICES

## ELASTIC MANDIBULAR ADVANCEMENT EMA DEVICE

Invented by Dr. Donald Frantz in 1993, the EMA is FDA cleared for the treatment of both obstructive sleep apnea and snoring. The EMA oral appliance from Myerson increases airway space by advancing the mandible using interchangeable elastic straps. EMA elastic come in 9 difference lengths, and 4 different strengths. The shorter the strap, the farther the mandible is advanced.

### EMA DESIGN ADVANTAGES INCLUDE:

- Metal free design.
- Easy to titrate in 1mm increments with EMA elastic straps.
- EMA elastic straps are latex free and available in 4 strengths (Soft, Medium, Firm and Extra Firm).



\$100  
(each)

## TAP DREAMTAP

The Thornton Adjustable Positioner® (TAP®) is a simple, patient-friendly treatment solution for snoring and sleep apnea. The TAP holds the lower jaw in a forward position so that it does not fall open during the night and cause the airway to collapse. It maintains a clear airway to reduce snoring and improve breathing. The unique design allows the

patient to fine-tune his/her treatment at home and work with the clinician to achieve the best results. Patients are empowered to manage the degree of lower jaw protrusion over as many nights as it takes to achieve the optimal treatment position. The TAP has a single point of central adjustment preventing an irregular bite and discomfort that may result from uneven bilateral adjustment.



\$315  
(each)

## SOMNODENT MAS CLASSIC ORIGINAL INNOVATION

- Highest quality acrylic
- Medical grade stainless steel ball clasps for retention
- Lingual-less option available for additional tongue space



\$300  
(each)

## SOMNODENT MAS HERBST ADVANCED VISUAL SIMPLICITY

- E0486 Medicare verified
- Visual calibration indicator
- 8mm range of advancement
- Reduces need for bite reset



\$300  
(each)

## SLIM-LINE OASYS ORAL/NASAL AIRWAY SYSTEM

Oral/nasal airway system is the first system to address the resistance to breathing in the entire upper airway. It is FDA approved as a nasal dilator

and mandibular repositioner for the treatment of both nasal and pharyngeal restrictions that cause snoring and sleep apnea.



\$400  
(each)

## PANTHERA DIGITAL SLEEP APNEA DEVICE

The next generation of digital sleep apnea devices that utilizes CAD/CAM technology to provide patients with:

- Smallest, least invasive and strongest FDA device on the market.
- Flexible & resistant to severe bruxing.
- Easily titrated with 1mm and .5mm advancement rods.
- Fabricated to HIP mounting and phonetic bite protocols.



\$175  
(each)

# DIGITAL/PRINTED SLEEP DEVICES



## DDSO POSTERIOR CONTACT

- Standard posterior pads
- Designed Digitally and 3D printed in white biocompatible Nylon
- Modular tongue buttons to promote the tongue to move up and forward

\$175 (each)

## DDSO WITH VERTICAL SHIMS

\$200 (each)

### VERTICAL SHIMS

- 1mm
- 2mm
- 3mm
- 4mm



# TMD DEVICES

## OLMOS DAY (OD) DECOMPRESSION

### PMT Thermoform

Indicated for chronic disc-displacement of skeletal and muscular asymmetry. Functions as positioning stent to correct vertical dimension, rotation, cant and protrusive changes. Lingual Anatomical ridge guides disclusion and protects contralateral joint in working movements.



\$125  
(each)

## NYLON ADJUSTMENT BURS APPLIANCE ADJUSTMENT

- Plain toothing with transverse section
- Used to relieve/remove material from Nylon devices



\$90  
(3 pack)

## OLMOS NIGHT DEPROGRAMMER

### Nylon Type12

- Control Parafunctional activity.
- Reduces the forces of clenching and grinding from chronic pain conditions.
- No suspected apnea.
- Upper and lower anterior pad so crowded or worn dentition are not a factor.
- No posterior contact.
- Full coverage to minimize tooth movement.
- Pressure formed (PMT).



\$150  
(each)

## OLMOS NIGHT POSITIONER

### Nylon Type12


- Combines benefits of Anterior pad with no posterior occlusion and ring prevents retrusion of the mandible in the supine position.
- Reduces inflammation (capsulitis).
- Prevents locking in supine position.
- Improve airway obstruction.
- Reduces forces of clenching or grinding.



\$165  
(each)

## UNIVERSAL POLISHING APPLIANCE FINISHING

- Scotch-Brite brushes for finishing and polishing acrylic or nylon devices
- Green: Fine Grit
- Black: Medium Grit
- Brown: Coarse Grit



\$20  
(10 pack)