

SOMNODENT SLEEP APNEA APPLIANCE

SOMNOMED®

offers a broad range of oral devices for the first line treatment of mild to moderate obstructive sleep apnea.

SOMNODENT DEVICES ARE:

COMFORTABLE FOR PATIENTS

- 88% of patients reported regular use of their SomnoDent device
- Unique SMH BFlex technology provides unparalleled comfort
- Patients can talk and open/close their mouths and lips while wearing device

HIGH QUALITY

- FDA 510K cleared and manufactured at an ISO 13485 certified facility
- Customized using the highest quality acrylic which does not discolor or attract odors
- Device life expectancy 3-5 years; and devices have a 3 year warranty against manufacturing defects

UNIQUELY COMPLIANT

- The 1st and only dorsal fin oral device with micro-recorder that objectively records compliance data
- Reports include hours device is worn and patient's supine or non-supine head position

CLINICALLY EFFECTIVE

- 15+ independent studies demonstrate the device's significant clinical benefits
- 91% of patients reported improvement in sleep quality while wearing SomnoDent

MEDICARE APPROVED

PDAC (Pricing, Data Analysis and Coding) contractor-verified for Medicare reimbursement:

Code: E0486

SOFT-LINER OPTIONS

- Comfortable for all-night wear
- Custom-fabricated for precise fit
- Provides a wide range of adjustment

COMFORT

QUALITY

COMPLIANT

EFFECTIVE

FUSION™

CUSTOM CALIBRATION

- Patient specific calibration
- 3 removable wings for 1mm of protrusion on the bottom
- Fine calibration in .1mm increments with lugs on the top
- Longer range of advancement (8.5mm) leading to fewer device resets



FLEX

PREMIUM COMFORT

- Proprietary SMH BFlex inner liner for superior retention
- Drop in fit - even for patients with short teeth, little undercut or crown / bridge work



CLASSIC

ORIGINAL INNOVATION

- Highest quality acrylic
- Medical grade stainless steel ball clasps for retention
- Lingual-less option available for additional tongue space



G2

LATEST INNOVATION

- Ease of calibration
- Wide range of calibration (-2mm to +8mm)
- Standard metal-free, also available with ball clasps



HERBST ADVANCE™

VISUAL SIMPLICITY

- E0486 Medicare verified
- Visual calibration indicator makes it easy to measure how far device has been advanced
- SMH BFlex offers superior comfort and retention
- Longer range of advancement (8mm) leading to fewer device resets over the life of the device



HERBST®

POPULAR VALUE

- E0486 Medicare verified
- Standard telescopic mechanism (also available in Shim)
- SMH BFlex offers superior comfort and retention, also available with ball clasps



SLEEP/ OSA APPLIANCES



TAP APPLIANCES

The Thornton Adjustable Positioner® (TAP®) is a simple, patient-friendly treatment solution for snoring and sleep apnea. The TAP holds the lower jaw in a forward position so that it does not fall open during the night and cause the airway to collapse. It maintains a clear airway to reduce snoring and improve breathing.

The unique design allows the patient to fine-tune his/her treatment at home and work with the clinician to achieve the best results. Patients are empowered to manage the degree of lower jaw protrusion over as many nights as it takes to achieve the optimal treatment position. The TAP has a single point of central adjustment preventing an irregular bite and discomfort that may result from uneven bilateral adjustment.



SnoreHookSplint SNOREHOOK APPLIANCE

